INSTRUCTION MANUAL

Electrophoretic Mini-Blotting Systems Double & Triple Wide

EBU-302

EBU-402


TABLE OF CONTENTS

	Page
Important User Information	3-4
Section 1 General Information	
1.1 Introduction	5
1.2 Specifications	5
1.3 Safety	5
Section 2 Description of parts	
2.1 Unpacking	6
2.2 Components/Assembly	6
Section 3 Instructions for Use	
3.1 Blotting Unit Preparation	7
3.2 Preparing the Gel for Transfer	7-8
3.3 Standard Electro-Blot Transfer	9
3.4 Removing the Cassettes	9
Section 4 Running Conditions	
4.1 Recommended Buffers & Running Conditions	10
4.2 References	10
Section 5 Maintenance of Equipment	
5.1 Care and Handling	11
5.2 Maintenance	11
Section 6 Equipment and Accessories	12

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first, and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. Scientific warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. Scientific warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. Scientific will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. Scientific shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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SECTION 1 General Information

1.1 Introduction

Electroblotting is a fast, easy, and efficient technique for transferring proteins from PAGE gels nitrocellulose and nucleic acids from agarose/PAGE gels onto nylon transfer membranes. The Wide Mini-Blotting Systems can effectively complete transfer in less than 2 hours compared to 6 to 24 for the more traditional "capillary" transfer. The C.B.S. Scientific Electrophoretic Blotting Systems come equipped with two electrode panels, one color-coded gel cassette assembly, holding tank with cooling chamber and exit ports and safety cover with attached power leads. The EBU-302 was designed primarily in conjunction with the Triple Wide Mini-Vertical Gel System (MGV-102/202-33), although the gel cassette will accommodate gels smaller than the set dimensions of 9cm (H) x 32.8cm (W) and can transfer up to two gels at one time. The EBU-402 was designed for use with the MGV-202-20 Double-Wide Mini-Vertical System, and can transfer two double wide gels with 20 x 10cm dimensions, or four standard mini gels.

1.2 Specifications

Constructions:

Buffer chamber, safety cover	Acrylic
Electrode panels	Acrylic and platinum wire .012" diameter
Power cords	FR Urethane rated 7500VDC @ 200mA, 65°C
Gel cassette	Perforated fiberglass panels red and black, two Scotch-Brite® pads, one sponge pad
Safety Certification	EN61010-1-1993 (IEC1010-1)

	EBU-402	EBU-302
Shipping Weight	7lbs	9lbs
Overall Size (l)x(w)x(h) cm	24x9.5x16	38x10x16
Maximum Gel Size - cm	20x9	32.8(w)x9(h)
Distance between electrodes -cm	5.3cm	5.5cm
Recommended buffer volume		
1 cassette	1.5L	2L
2 cassettes	1.4L	1.9L

1.3 Safety

Power to the Electro-Blotter is to be supplied by an external DC voltage power supply that must be ground isolated so that the DC voltage output floats with respect to ground. For any power supply used, the maximum specified operating parameters for the units are:

Maximum Limits

100 VDC
150 watts power
1500mA current
50°C ambient temperature

Current to the unit, provided from the external power supply, must enter the unit through the lid assembly, providing a safety interlock to the user. Current to the unit is broken when the lid is removed. **Do not attempt to use the unit without the safety lid, and always turn the power supply off before removing the lid, or when working with the unit in any way. Follow safety precautions specified by the power supply manufacturer.**

SECTION 2

Description of Parts


2.1 Unpacking the EBU-302 or EBU-402

Please verify that your unit comes complete with the following components:

Blotting chamber

- Blotting cassette with sponge pads
- Safety cover with attached power leads
- 2 electrode panels
- Tubing adapters for cooling base

2.2 Components/Assembly

SECTION 3

Instructions for Use

3.1 Blotting Unit Preparation

1. Place the blotting chamber on a level work surface in an authorized work area.
2. Consult applications table 4.1 for type of transfer, buffer system, membrane type and power settings.

Coolant Circulation:

For coolant circulation, connect a heavy wall tubing which will not kink (hose clamps are optional), to the tubing adapters supplied with the unit (excessive pressure will damage the seal between the glass and the base and possibly create leaks). Attach a regulated pump or recirculating water bath (follow manufacturer's instructions), recommended flow is 500 ml per minute or 15psi. **Do NOT** use tap or house water as it can be subject to large fluctuations in pressure.

1. Place, do NOT drop, a magnetic spin bar (not supplied by C.B.S.) in the bottom of the tank. This maintains pH, buffer and heat circulation during the electroblot.
2. Place tank on top of the magnetic stirrer.
3. Connect to a cold water supply (see coolant circulation above).
4. Buffer should be pre-cooled to 10°C.

3.2 Preparing the Gel for Transfer

The EBU-302 and EBU-402 can hold up to 2 gel cassettes simultaneously with space provided for a magnetic stir bar to allow buffer circulation and heat exchange for uniform transfers. When the resolving gel electrophoresis is complete, proceed with staining and photo-documentation if applicable. Cut one corner off the gel so that correct orientation is maintained throughout the procedure.


Preparation of Electroblotting Components

1. Cut transfer membrane to size of gel or sufficient dimensions to cover the relevant bands. Also, pre-cut Whatman-type 3MM filter paper to same size and soak in transfer buffer until completely saturated (15 to 30 minutes).
2. Pre-wet membrane chosen according to manufacturer's instructions. Nylon or Nylon-supported nitrocellulose should be soaked in ddH₂O. PVDF in MeOH. Equilibrate all types in transfer buffer for 3 minutes.
3. Pre-equilibrate gel in transfer buffer to be used prior to electroblotting for 10 to 30 minutes depending on gel thickness.
4. Open the transfer cassette and submerge the red (+) panel in a shallow plastic or glass (Pyrex) dish. Fill with enough buffer to cover entire cassette. Submerge a single Sponge® pad and gently tease to release trapped air bubbles until pad is completely saturated.
5. Layer on top of the Sponge pad your transfer membrane. (Some prefer to layer 1 to 2 sheets of saturated 3mm filter paper). Carefully search for and remove trapped air bubbles at each layer by rolling a pipet over the surface.


Note: Sponge pads will degrade/compress over time. To maintain tightness in the stack, add filter paper or replace pads.
6. Apply gel correctly oriented to submerged transfer membrane. Check for air pockets and remove. Note: Some proteins will begin transfer immediately upon contracting the transfer membrane. Disturbing or moving the gel/membrane interface can result in a smeared blot.

7. Quickly follow with 1-2 sheets of saturated 3mm filter paper (optional).
8. Place second saturated Sponge pad (referenced below as Scotch Brite) over sandwich assembly and close the cassette.

Cassette Assembly

9. Remove entire assembly from pyrex dish and load into desired position (slot) in buffer chamber with black screen facing towards black terminal (cathode) and red screen facing towards red terminal (anode).

3.3 Standard Electro-Blot Transfer

1. Fill the chamber with transfer buffer. (Buffers should be prepared fresh with reagent grade chemicals and pre-cooled). The transfer buffer should be pre-cooled and added to the top of the platinum labyrinth. The buffer should not come in contact with the banana plugs when the gel cassette sandwich is immersed in the unit.
2. Align safety cover over the unit and carefully attach. Begin cycling of coolant.
3. Connect the leads to the power supply, matching the color-coded red to red and black to black. **See Section 4.1 for recommended power conditions.**
4. Transfer times will vary according to several parameters. Optimization of electro-blotting transfers must be determined empirically. Keep in mind the following principles that govern the movement of molecules of gel electrophoresis:

-Thicker or higher percentage gels will take longer to transfer than thinner or lower percentage gels.

-Large molecules will need extended transfer times to completely transfer.

-Actual transfer times for defined conditions can be approximated by running molecular weight standards.

3.4 Removing the Cassettes

1. Turn the power supply off and disconnect the leads from the power supply.
2. Remove the safety cover from the unit, by placing thumbs on white posts next to red & black connectors, then pushing down while pulling up with fingers under lid. **DO NOT pull on power cords.**
3. Gently lift the cassette from the unit. **Always wear gloves, eye protection and protective clothing** if buffer and/or gel contain Ethidium Bromide. Ethidium Bromide is a powerful mutagen; gloves, eye protection and protective clothing should always be worn when handling the gel or buffer solutions. See Material Data Safety Sheets if using various dye solutions.
4. Mark the orientation of the membrane with a pencil or by cutting off a corner and take apart the cassette carefully.
5. Process membrane according to type of transfer and manufacturer's recommendations

SECTION 4

Applications & Running Conditions for Tank Type Electro-Blotting

4.1 Recommended Buffers, Power Settings, and Transfer Times

Precise transfer conditions will vary according to the number and type of gels used, buffer conditions employed, power input, and the general goal of the experiment. Refer to the reference section for in depth discussions on practical and theoretical approaches to gel transfers.

Using the enclosed Towbin buffer system use the following conditions at constant voltage:

At 4°C: 26V, current range: 250-300mA for 18 hours

At room temperature: 26V, current range: 400-500mA for 18 hours

Towbin Buffer – 1X, 20% MeOH:

0.025M Tris Base

0.192 M Glycine

0.05-0.1%(w/v) SDS

20% (v/v) Methanol

4.2 References

1. Ausubel, F.M., Brent, R., Kingston, R.E., Moore, D.D., Seidman, J.G., Smith, J.A., Struhl, K. (ed.) (1993). *Current Protocols in Molecular Biology*. Vol. 2, Greene Publishing Associates, Inc. and John Wiley & Sons, Inc., Ch.10.
2. Burnette, W.N. (1981). Western blotting: Electrophoretic transfer of proteins from sodium dodecyl sulfate-polyacrylamide gels to unmodified nitrocellulose and radiographic detection with antibody and radioiodinated protein A. *Anal. Biochem.* 112:195-203.
3. Peluso, R.W. and Rosenberg, G.H. (1987). Quantitative electrotransfer of proteins from sodium dodecyl sulfate polyacrylamide gels onto positively charged nylon membranes. *Anal. Biochem.* 162:389-398.
4. Perides, G. Plagens, U., and Traub, P. (1986). Protein transfer from fixed, stained and dried polyacrylamide gels and immunoblot with protein A-gold. *Anal. Biochem.* 152:94-99.
5. Tesfaigzi, J., Smith-Harrison, W., and Carlson, D.M. (1994). A simple method for reusing western blots on PVDF membranes. *Biotechniques.* 17:268-269.
6. Towbin, H., Staehelin, T., and Gordon, J. (1979). Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: Procedure and some applications. *Proc. Natl. Acad. Sci. U.S.A.* 76: 4350-4354.
7. Sambrook, J., Fritsch, E.F., Maniatis, T. (1989). *Molecular Cloning. A Laboratory Manual*. 2nd ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York. 18.47-18.61.

SECTION 5

Maintenance of Equipment

5.1 Care and Handling

The plastic components of the Electrophoretic Blotting units are fabricated from acrylic and polycarbonate. Electrodes and connectors are made from pure platinum, stainless steel, and chrome plated brass. As with any laboratory instrument, adequate care ensures consistent and reliable performance.

After each use, rinse buffer chamber, gel tray and combs with de-ionized water. Wipe dry with a soft cloth or paper towel, or allow to air dry. Whenever necessary, all components may be washed gently with water and a non-abrasive detergent, and rinsed and dried as above. *Never* use abrasive cleaners, glass cleaning sprays or scouring pads to clean the components, as these will damage the unit and components.

Additional precautions:

- Do not autoclave or dry-heat sterilize the apparatus or components.
- Do not expose the apparatus or components to phenol, acetone, benzene, halogenated hydrocarbon solvents or alcohols.
- Avoid prolonged exposure of the apparatus or components to UV light.
- For Northern blotting, DO NOT treat with diethylpyrocarbonate (DEPC)-treated water for extended periods at 37°C.
- Electrostatic charge and heat affix the conductive platinum surface covering the titanium. It can be scratched off or damaged by using sharp objects or scouring pads. This will disturb the continuity of the electric field and have a negative effect on transfers.
- The stainless electrode can be permanently damaged by corrosion if the polarity between anode and cathode is reversed. CAUTION: Be certain of power connections before initiating transfer.

5.2 Maintenance

The following inspection and maintenance procedures will help maintain the safety and reliable performance of the Blotting systems. Replacement parts can be ordered by calling 1-858-755-4959 or by contacting your local distributor.

- Banana plugs and power cords should be inspected regularly. If the banana plugs become loose or do not feel friction tight replace the plugs or power cords.
- Should power cord assemblies (connectors, wire or shrouds) show any signs of wear or damage (e.g. cracks, nicks, abrasions, or melted insulation), replace them immediately.
- The platinum wire is secured to the banana jack by compression between a stainless washer and the jack nut. The nut/washer interface should be tight and free of corrosion.

CONTACT US

Telephone:
Local or International
858-755-4959
Toll Free: 858-755-0733

Fax: 858-755-0733

Check out our new and improved online ordering!

www.cbsscientific.com

E-mail address:
sales@cbssci.com

Mailing address:
C.B.S. Scientific Company
P.O. Box 856
Del Mar, CA 92014

Shipping address:
C.B.S. Scientific Company
10805 Vista Sorrento Pkwy
Suite 100
Solana Beach, CA 92121

Credit Card Options:
Visa/Mastercard
Discover/American Express