INSTRUCTION MANUAL

Electrophoretic Blotting System

EBU-102
EBU-102P


TABLE OF CONTENTS

	Page
Important User Information	3-4
Section 1 General Information	
1.1 Introduction	5
1.2 Specifications	5
1.3 Safety	5
Section 2 Description of parts	
2.1 Unpacking EBU-102 or EBU-102P	6
2.2 Components/Assembly	6
2.3 Unpacking EBU-202	7
2.4 Components/Assembly	7
Section 3 Instructions for Use	
3.1 Blotting Unit Preparation	8
3.2 Preparing the Gel for Transfer	8-9
3.3 Standard Electro-Blot Transfer	10
3.4 Removing the Cassettes	10
Section 4 Running Conditions	
4.1 Recommended Buffers, Membranes & Running Conditions	11
4.2 References	12
Section 5 Maintenance of Equipment	
5.1 Care and Handling	13
5.2 Maintenance	13
Section 6 Equipment and Accessories	14

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first, and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. Scientific warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. Scientific warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. Scientific will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. Scientific shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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SECTION 1 General Information

1.1 Introduction

Electroblotting is a faster, easier, and more efficient method to permanently transfer PAGE gels to nitrocellulose and nucleic acids to nylon, than capillary (passive diffusion) blotting. For a similar transfer, traditional "capillary blotting" can take 6-24 hours to complete compared to only 2 hours or less using electroblotting. The C.B.S. Scientific Electrophoretic Blotting Systems come equipped with; two electrode panels, one color-coded gel cassette assembly, electrophoresis tank with cooling chamber and safety cover with attached power leads. The EBU-102 or EBU-102P gel cassette size is 15cm x 21.5cm and can transfer up to three gels at one time. The EBU-202 gel cassette size is 9 x 11cm and can transfer up to 4 gels at one time.

1.2 Specifications

Constructions:

Buffer chamber, safety cover	Acrylic
Electrode panels	Acrylic and platinum wire .012" diameter
Power cords	FR Urethane rated 7500VDC @ 200mA, 65°C
Gel cassette	Perforated fiberglass panels red and black, two Scotch-Brite® pads, and one sponge pad
Safety Certification	EN61010-1-1993 (IEC1010-1)

	EBU-102 or EBU-102P	EBU-202
Shipping Weight	7 lbs.	5 lbs.
Overall Size (l)x(w)x(h) cm	19 x 11 x 28.5	14 x 13 x 16
Maximum Gel Size - cm	15.3 x 21.5	10 x 10
Distance between electrodes -cm	8.7	10
Recommended buffer volume		
1 cassette	3.75 liters	1.1 liter
2 cassettes	3.76 liters	0.9 liter
3 cassettes	3.77 liters	0.8 liter
4 cassettes	N/A	0.7 liter

1.3 Safety

Power to the Electro-Blotter is to be supplied by an external DC voltage power supply that must be ground isolated so that the DC voltage output floats with respect to ground. For any power supply used, the maximum specified operating parameters for the units are:

Maximum Limits

100 VDC
150 watts power
1500mA current
50°C ambient temperature

Current to the unit, provided from the external power supply, must enter the unit through the lid assembly, providing a safety interlock to the user. Current to the unit is broken when the lid is removed. **Do not attempt to use the unit without the safety lid, and always turn the power supply off before removing the lid, or when working with the unit in any way. Follow safety precautions specified by the power supply manufacturer.**

SECTION 2


Description of Parts

2.1 Unpacking the EBU-102 or EBU-102P

Please verify that your unit comes complete with the following components:

- Blotting chamber
- Blotting cassette with sponge pads
- Safety cover with attached power leads
- 2 electrode panels
- Tubing adapters for cooling base

2.2 Components/Assembly


2.3 Unpacking the EBU-202

Please verify that your unit comes complete with the following components:

- Blotting chamber
- Blotting cassette with Scotch-Brite pads
- Safety cover with attached power leads
- 2 electrode panels
- Tubing adapters for cooling base

2.4 Components/Assembly

SECTION 3

Instructions for Use

3.1 Blotting Unit Preparation

1. Place the blotting chamber on a level work surface in an authorized work area.
2. Consult applications table 4.1 for type of transfer, buffer system, membrane type and power settings.

Coolant Circulation:

For coolant circulation, connect a heavy wall tubing which will not kink (hose clamps are optional), to the tubing adapters supplied with the unit (excessive pressure will damage the seal between the glass and the base and possibly create leaks). Attach a regulated pump or recirculating water bath (follow manufacturer's instructions), recommended flow is 500 ml per minute or 15psi. **Do NOT** use tap or house water as it can be subject to large fluctuations in pressure.

1. Place, do NOT drop, a magnetic spin bar (not supplied by C.B.S.) in the bottom of the tank. This maintains pH, buffer and heat circulation during the electroblot.
2. Place tank on top of the magnetic stirrer.
3. Connect to a cold water supply (see coolant circulation above).
4. Buffer should be pre-cooled to 10°C.

3.2 Preparing the Gel for Transfer

The EBU-202 can hold up to 4 gel cassettes simultaneously with space provided for a magnetic stir bar to allow buffer circulation and heat exchange for uniform transfers. The EBU-102 holds up to 3 gel cassettes simultaneously. When the resolving gel electrophoresis is complete, proceed with staining and photo-documentation if applicable. Cut one corner off the gel so that correct orientation is maintained throughout the procedure.

Preparation of Electroblotting Components

1. Cut transfer membrane to size of gel or sufficient dimensions to cover the relevant bands. Also, pre-cut Whatman-type 3MM filter paper to same size and soak in transfer buffer until completely saturated (15 to 30 minutes).
2. Pre-wet membrane chosen according to manufacturer's instructions. Nylon or Nylon-supported nitrocellulose should be soaked in ddH₂O. PVDF in MeOH. Equilibrate all types in transfer buffer for 3 minutes.
3. Pre-equilibrate gel in transfer buffer to be used prior to electroblotting for 10 to 30 minutes depending on gel thickness.
4. Open the transfer cassette and submerge the red (+) panel in a shallow plastic or glass (Pyrex) dish. Fill with enough buffer to cover entire cassette. Submerge a single Scotch Brite® pad and gently tease to release trapped air bubbles until pad is completely saturated.
5. Layer on top of the Scotch Brite pad your transfer membrane. (Some prefer to layer 1 to 2 sheets of saturated 3mm filter paper). Carefully search for and remove trapped air bubbles at each layer by rolling a pipet over the surface.


Note: Scotch Brite pads will degrade/compress over time. To maintain tightness in the stack, add filter paper or replace pads.

6. Apply gel correctly oriented to submerged transfer membrane. Check for air pockets and remove. Note: Some proteins will begin transfer immediately upon contracting


the transfer membrane. Disturbing or moving the gel/membrane interface can result in a smeared blot.

7. Quickly follow with 1-2 sheets of saturated 3mm filter paper (optional).
8. Place second saturated Scotch-Brite pad over sandwich assembly and close the cassette.

Cassette Assembly

9. Remove entire assembly from pyrex dish and load into desired position (slot) in buffer chamber with black screen facing towards black terminal (cathode) and red screen facing towards red terminal (anode).

3.3 Standard Electro-Blot Transfer

1. Fill the chamber with transfer buffer. (Buffers should be prepared fresh with reagent grade chemicals and pre-cooled). The transfer buffer should come up to the top of the platinum labyrinth. The buffer should not come in contact with the banana plugs when the gel cassette sandwich(es) are immersed in the unit.
2. Align safety cover over the unit and carefully attach. Begin cycling of coolant.
3. Connect the leads to the power supply, matching the color-coded red to red and black to black. **See Section 4.1 for recommended power conditions.**
4. Transfer times will vary according to several parameters. Optimization of electro-blotting transfers must be determined empirically. Keep in mind the following principles that govern the movement of molecules of gel electrophoresis:
 - Thicker or higher percentage gels will take longer to transfer than thinner or lower percentage gels.
 - Large molecules will need extended transfer times to completely transfer.
 - Actual transfer times for defined conditions can be approximated by running molecular weight standards.

3.4 Removing the Cassettes

1. Turn the power supply off and disconnect the leads from the power supply.
2. Remove the safety cover from the unit, by placing thumbs on white posts next to red & black connectors, then pushing down while pulling up with fingers under lid. **DO NOT pull on power cords.**
3. Gently lift the cassette from the unit. **Always wear gloves, eye protection and protective clothing** if buffer and/or gel contain Ethidium Bromide. Ethidium Bromide is a powerful mutagen; gloves, eye protection and protective clothing should always be worn when handling the gel or buffer solutions. See Material Data Safety Sheets.
4. Mark the orientation of the membrane with a pencil or by cutting off a corner and take apart the cassette carefully.
5. Process membrane according to type of transfer and manufacturer's recommendations.

SECTION 4

Applications & Running Conditions for Tank Type Electro-Blotting

4.1 Recommended Buffers, Membranes, Power Settings and Transfer Times

Gel Type	Blotting Buffer	Transfer Membrane	Power
denaturing: SDS-PAGE (Western)	Towbin (25mM Tris/192mM Glycine/ pH 8.3/Methanol 20% (w/v)	nitrocellulose/nylon PVDF: ion-exchange	<ul style="list-style-type: none"> room temperature, 70V-100V, 200mA, 1.5- 4 hours 4°C, 150V, 400mA, .5-2 hours Overnight, 30V, 75mA
non-denaturing: SDS-PAGE (acidic or neutral proteins)	25mMTris/192mM Glycine, pH 8.3 25mM sodium phosphate, pH6.5 15mM sodium borate, pH9.2	nitrocellulose; nylon PVDF; ion-exchange diazo-paper diazo-paper	<ul style="list-style-type: none"> room temperature, 60V, 200mA, 4 hours 4°C, 125V, 400mA, 1.5 hours Overnight, 30V, 100mA
IEF (native- basic proteins, acid urea gels)	0.7% acetic acid (w/v)*	nitrocellulose; nylon diazo-paper	<ul style="list-style-type: none"> room temperature, 70V, 400mA, 4 hours 4°C, 100V, 600mA, 1.5 hours Overnight, 30V, 200mA
DNA (Southern)	1XTAE or agarose 0.5X TBE, acrylamide	nylon; +charged nylon, neutral or +charged	<ul style="list-style-type: none"> room temperature, 30-40V, 25-150mA, 2-4 hours 4°C, 125V, 400mA, 1.5 hours
RNA (Northern)	19.6mM phosphate/5.4mM citrate pH3.0 or 25mM sodium phosphate pH 6.5	nylon nitrocellulose	<ul style="list-style-type: none"> room temperature, 60V, 200mA, 4 hours 4°C, 125V, 400mA, 1.5 hours Overnight, 30V, 200mA

4.2 References

- Ausubel, F.M., Brent, R., Kingston, R.E., Moore, D.D., Seidman, J.G., Smith, J.A., Struhl, K. (ed.) (1993). *Current Protocols in Molecular Biology*. Vol. 2, Greene Publishing Associates, Inc. and John Wiley & Sons, Inc., Ch.10.
- Burnette, W.N. (1981). Western blotting: Electrophoretic transfer of proteins from sodium dodecyl sulfate-polyacrylamide gels to unmodified nitrocellulose and radiographic detection with antibody and radioiodinated protein A. *Anal. Biochem.* 112:195-203.
- Peluso, R.W. and Rosenberg, G.H. (1987). Quantitative electrotransfer of proteins from sodium dodecyl sulfate polyacrylamide gels onto positively charged nylon membranes. *Anal. Biochem.* 162:389-398.
- Perides, G. Plagens, U., and Traub, P. (1986). Protein transfer from fixed, stained and dried polyacrylamide gels and immunoblot with protein A-gold. *Anal. Biochem.* 152:94-99.
- Tesfaigzi, J., Smith-Harrison, W., and Carlson, D.M. (1994). A simple method for reusing western blots on PVDF membranes. *Biotechniques.* 17:268-269.
- Towbin, H., Staehelin, T., and Gordon, J. (1979). Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: Procedure and some applications. *Proc. Natl. Acad. Sci. U.S.A.* 76: 4350-4354.
- Sambrook, J., Fritsch, E.F., Maniatis, T. (1989). *Molecular Cloning. A Laboratory Manual*. 2nd ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York. 18.47-18.61.

SECTION 5 Maintenance of Equipment

5.1 Care and Handling

The plastic components of the Electrophoretic Blotting units are fabricated from acrylic and polycarbonate. Electrodes and connectors are made from pure platinum, stainless steel, and chrome plated brass. As with any laboratory instrument, adequate care ensures consistent and reliable performance.

After each use, rinse buffer chamber, cassettes and Scotch Brite pads with de-ionized water. Wipe dry with a soft cloth or paper towel, or allow to air dry. Whenever necessary, all components may be washed gently with water and a non-abrasive detergent, and rinsed and dried as above. *Never* use abrasive cleaners, glass cleaning sprays or scouring pads to clean the components, as these will damage the unit and components.

Additional precautions:

- Do not autoclave or dry-heat sterilize the apparatus or components.
- Do not expose the apparatus or components to phenol, acetone, benzene, halogenated hydrocarbon solvents or alcohols.
- Avoid prolonged exposure of the apparatus or components to UV light.
- Do NOT treat with diethylpyrocarbonate (DEPC)-treated water for extended periods at 37°C. A brief rinse with DEPC-water is sufficient after a thorough wash, followed by a quick rinse in 50% ethanol.
- Electrostatic charge and heat affix the conductive platinum surface covering the titanium. It can be scratched off or damaged by using sharp objects or scouring pads. This will disturb the continuity of the electric field and have a negative effect on transfers.
- The stainless electrode can be permanently damaged by corrosion if the polarity between anode and cathode is reversed. CAUTION: Be certain of power connections before initiating transfer.

5.2 Maintenance

The following inspection and maintenance procedures will help maintain the safety and reliable performance of the Blotting systems. Replacement parts can be ordered by calling 1-858-755-4959 or by contacting your local distributor.

- Banana plugs and power cords should be inspected regularly. If the banana plugs become loose or do not feel friction tight replace the plugs or power cords.
- Should power cord assemblies (connectors, wire or shrouds) show any signs of wear or damage (e.g. cracks, nicks, abrasions, or melted insulation), replace them immediately.
- The platinum wire is secured to the banana jack by compression between a stainless washer and the jack nut. The nut/washer interface should be tight and free of corrosion.

SECTION 6

Equipment and Accessories

Cat. #	Item
EBU-102	Electrophoretic Blotting System. Includes chamber with 2 platinum wire electrode panels, 1 color-coded gel cassette, safety cover with power leads
EBU-102P	Electrophoretic Blotting System, planar electrodes. Includes chamber with 2 planar electrode panels, 1 color-coded gel cassette, safety cover with power leads electrode panels, safety cover and interlocking power leads.
EBE-102	Replacement Electrode Panel for EBU-102
EBE-102P	Replacement Electrode Panel for EBU-102P
EBE-202	Replacement Electrode Panel for EBU-202
EBC-102	Additional Gel Cassette, color coded
EBS-102	Sponge Pads

POWER SUPPLIES

Cat. #	Item
EPS-200-II	High Current Power Supply with timer, 110/50-60Hz CE. CV or CC 5-200V, current range: 4-1500mA, 200 Watts. -200 Volts, 4-500mA, 90 Watts, with timer, 110V/60Hz
EPS-200-IIV	High Current Power Supply with timer, 220/50Hz CE. CV or CC 5-200V, current range: 4-1500mA, 200 Watts. -200 Volts, 4-500mA, 90 Watts, with timer, 110V/60Hz
EPS-300-II	Mini Power Supply, 10-300 Volts, 4-500mA, 90 Watts, with timer, 110V/60Hz
EPS-300-IIV	Mini Power Supply, 10-300 Volts, 4-500mA, 90 Watts, with timer, 220V/50Hz

CONTACT US

Telephone:
Local or International
858-755-4959
Toll Free: 858-755-0733

Fax: 858-755-0733

Check out our new and improved online ordering!

www.cbsscientific.com

E-mail address:
sales@cbssci.com

Mailing address:
C.B.S. Scientific Company
P.O. Box 856
Del Mar, CA 92014

Shipping address:
C.B.S. Scientific Company
10805 Vista Sorrento Pkwy
Suite 100
Solana Beach, CA 92121

Credit Card Options:
Visa/Mastercard
Discover/American Express