INSTRUCTION MANUAL

Horizontal Midi-Gel Productivity Package

GCMGU-402T

GCMGU-502T

GCMGU-602T


TABLE OF CONTENTS

	Page
Important User Information	3-4
Section 1 General Information	
1.1 Introduction	5
1.2 Specifications	5
1.3 Safety	6
Section 2 Description of parts	
2.1 Unpacking	6
2.2 Components/Assembly	6
Section 3 Instructions for Use	
3.1 Midi-Gel Unit Preparation and Gel Casting	7-8
3.2 Running the Gel	8
3.3 Removing the Gel	9
3.4 Maximum Well/Comb volumes	10
Section 4 Running Conditions	
4.1 Recommended Power	10
4.2 Recommended Buffers	10-11
Section 5 Maintenance of Equipment	
5.1 Care and Handling	11
5.2 Maintenance	11
Section 6 Equipment and Accessories	12-13

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first, and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. Scientific warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. Scientific warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. Scientific will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. Scientific shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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General Information

1.1 Introduction

The mid-range Horizontal Midi-Gel units, MGU-402T (11 x 11cm), MGU-502T (10 x 14cm) and MGU-602T (14 x 14cm) effectively allow rapid separation and identification of various fragments from DNA, RNA, PCR™ products and synthetic oligonucleotides. When used in combination with two or more sample combs, designed for multi-channel pipette loading, multiple screening of large sample numbers can be facilitated easily. The unique tapered baffle system allows immediate gel casting directly onto removable trays without having to seal with agarose plugs or tape. The extensive variety of comb formats and material types ensure a fast and easy method of identifying critically important samples while using minimal quantities of expensive agarose.

Two methods for gel formation: Option 1, in Situ tapered baffle systems for simplified gel casting. The tapered baffle system allows gel casting directly onto the unit or removable trays without having to seal with agarose plugs or tape. Option 2: External horizontal gel casters. Precast your agarose gels for future use in an easy to use and easy to store fashion. These support trays have leveling bases and are lined with soft, coated neoprene to provide firm yet supple support for casting gels without leakage.

1.2 Specifications

Constructions:

Buffer chamber	Acrylic
Safety Cover	Acrylic
Electrodes	Platinum wire .012" diameter
Power cords	FR Urethane rated 7500VDC, 200mA, 65°C
Tapered baffles	Acrylic
Combs	Polycarbonate or Teflon
Adjustable comb backing	Acrylic
UVT Gel Tray	UV transparent acrylic
Safety Certification	EN61010-1-1993 (IEC1010-1)

	MGU-402T	MGU-502T	MGU-602T
Shipping weight	4lb.	6lb.	7lb.
Overall Size (l)x(w)x(h) cm	22.5 x 13 x 14	21.5 x 16.7 x 14	27.5 x 16.7 x 14
Gel Size cm	11 x 11	14 x 10	14 x 14
Recommended buffer volume: with gel tray	550 ml	700 ml	750 ml
Distance between electrodes	20 cm	19 cm	25 cm
Voltage limit	250VDC	250VDC	250VDC

1.3 Safety

Power to the MGU-402T, MGU-502T or MGU-602T unit is to be supplied by an external DC voltage power supply that must be ground isolated, so that the DC voltage output floats with respect to ground. For any power supply used, the maximum specified operating parameters for the units are:

Maximum Limits

250VDC voltage

15 watts power

60 mA current

55°C ambient temperature

Current to the unit, provided from the external power supply, must enter the unit through the lid assembly, providing a safety interlock to the user. Current to the unit is broken when the lid is removed. **Do not attempt to use the unit without the safety lid, and always turn the power supply off before removing the lid, or when working with the unit in any way. Follow safety precautions specified by the power supply manufacturer.**

SECTION 2


Description of Parts

2.1 Unpacking

Please verify that your unit comes complete with the following components:

- Midi-Horizontal electrophoresis unit
- Tapered baffles (2)
- Safety cover with attached DC power leads
- UV Transparent slotted gel tray (2)
- Two combs with adjustable height backing
- One gel casting tray, single place (not pictured)

2.2 Components/Assembly

SECTION 3 Instructions for Use

3.1 Midi-Gel Unit Preparation and Gel Casting



Place Midi-Gel unit on a level work surface in an authorized work area. Place the adjustable height comb assembly in the gel tray. Loosen the white screws holding the comb to the backing and adjust comb depth (generally 2mm). Gently tighten the adjusting screws and set the comb aside. **See Table 1 for comb volumes.**


The Midi-Gel unit provides the user with two options for gel casting: (a) **gel casting directly into unit using gel tray** or (b) **gel casting in gel tray outside of unit.**

A. Casting Directly into Unit Using Gel Tray.


1. Place gel tray onto raised bed of unit. Insert the tapered baffles, simultaneously on either side of gel tray and slide into the grooves next to the bed, handles to the outside. The tapered baffles will rest friction tight against the gel tray. **DO NOT** press baffles too tightly into stops.

2. Position the gel comb(s). Remembering, current travels from negative (black cathode) towards positive (red anode). Pour the agarose solution onto the gel bed. Allow the agarose slab to cool uncovered until it turns opaque. Agarose around combs will take a few more minutes to cool than other parts of the slab.


3. Slide the baffles out of their slots and store in the slots provided at either end of the unit near both terminals.


B. Gel Casting into Gel Tray Outside of Unit.

1. Tape both ends of gel tray with tape.

2. Place gel comb(s) into desired slot(s) in the gel tray. Using a pipette or a pipettor, run a bead of agarose solution along the edges of the gel tray where it meets the tape to assure a leak-proof seal. Then pour the agarose solution onto the gel tray. Allow the agarose slab to cool uncovered until it turns opaque. NOTE: (Optional storage) Pre-cast gels can be stored for up to one month if wrapped in plastic wrap and stored at 4°C. Remove combs prior to storing.

3. Remove gel tape at each end. Position the tray in the unit, remembering, current travels from negative (black cathode) towards positive (red anode).


C. Gel Casting with External Gel Caster.

1. Place external gel caster on level surface and level using thumb screws.
2. Push gel tray against one side of gel caster, then the other. Confirm it the gel tray is flat to the bottoms of the caster


3.2 Running the gel


1. Add enough buffer to fill both reservoirs and overflow the surface of the gel to a depth of 2-3mm. Gently remove comb(s). Flush out any air bubbles in the wells. Load the samples into the sample wells. NOTE: DO NOT FORGET TO LOAD DNA SIZE STANDARD.

2. Align safety cover over the unit and carefully attach, so not to disturb samples.


3. Connect the leads to the power supply, matching the color-coded red to red and black to black. **See Section 4.1 for recommended power conditions.** Begin separation by electrophoresis.

3.3 Removing the Gel

1. Turn the power supply off and disconnect the leads from the power supply.
2. Remove the safety cover from the unit, by placing thumbs on white posts next to red & black connectors, then pushing down while pulling up with fingers under lid. **DO NOT pull on power cords.**

3. Gently lift the gel tray from the unit. **Always wear gloves, eye protection and protective clothing** if buffer and/or gel contains Ethidium Bromide. Ethidium Bromide is a powerful mutagen, gloves, eye protection and protective clothing should always be worn when handling the gel or buffer solutions. See Material Data Safety Sheets.
4. View separated fragments under UV light, using proper protection for eyes and skin (see manufacturer's instructions).

3.4 Maximum Well/Comb Volumes

NOTE: To calculate sample well volume expressed in millimeters (mm) of height, divide maximum volume by tooth depth.

MBU Comb for MGU-402T # of wells	Tooth Depth mm	Overall length of comb mm	Spacing between teeth mm	Tooth width/mm	1.0mm thickness volume per tooth microliters (ul)	1.5mm thickness volume per tooth microliters (ul)	2.0mm thickness volume per tooth microliters (ul)	3.0mm thickness volume per tooth microliters (ul)
5	12.7	87.28	1.58	16.2	205.74	308.61	411.48	617.22
12	12.7	87.28	1.58	5.83	74.04	111.06	148.08	222.12
18	12.7	87.28	1.58	3.37	42.79	64.19	85.58	128.37
SGC14 Comb for MGU-502T & MGU-602T # of wells	Tooth Depth mm	Overall length of comb mm	Spacing between teeth mm	Tooth width/mm	1.0mm thickness volume per tooth microliters (ul)	1.5mm thickness volume per tooth microliters (ul)	2.0mm thickness volume per tooth microliters (ul)	3.0mm thickness volume per tooth microliters (ul)
1	12.7	120.64	3	120	1524	2286	3048	4672
3E	12.7	120.64	3	38.1	483.87	725.81	967.74	1451.61
5	12.7	120.64	3	21.6	274	411.48	548	822
10	12.7	120.64	3	9.2	117	175.26	232	350
12	12.7	120.64	2.4	7.8	99	148.59	198	297
13	12.7	120.64	2.4	6.6	83.82	125.73	167.64	251.46
16	12.7	120.64	2.4	5.3	67.3	100.9	134.6	201
20	12.7	120.64	1.6	4.5	57	85.7	114	171
26	12.7	120.64	1.6	3.1	39.3	59	78	118
27	12.7	120.64	1.6	2.93	37.21	55.82	74.42	111.63

SECTION 4

Running Conditions

4.1 Recommended Power

As a rule, optimal resolution of larger molecules is achieved during longer runs at lower voltages, whereas smaller molecules require shorter runs at higher voltages. Applied voltage gradients can therefore be anywhere in the range of 1-10 VDC/cm of gel. Using the standard buffer systems listed below, most runs will use 5 VDC/cm as a general rule.

The usual run time will vary for the voltage chosen but should range from 15 to 60 minutes. Nucleic acid migration is monitored by the progress of marker dyes. Constant power is not a necessity, but it produces uniform heat throughout the run, therefore minimizing band diffusion. Be sure the polarity is correct i.e. that the DNA is loaded near the cathode (black electrode) to run toward the anode (red terminal).

Agarose gels may be stored for several days at 4°C wrapped in plastic wrap. Seakem Agarose (FMC) is used (normally) for preparative and analytical gels. Other types of agarose can be used for special purposes.

4.2 Recommended Buffers*

Type*	Concentrated Stock/liter	Final Concentration
TAE (Tris-acetate)	50X - 242 gm Tris base 57.1ml glacial acetic acid 100 ml 0.5M EDTA (pH 8.0)	1X - 0.04M Tris-acetate 0.001M EDTA
TBE (Tris-borate)	5X - 54 gm Tris base 27.5 gm boric acid 20 ml 0.5M EDTA (pH 8.0)	0.5X - 0.045M Tris-borate 0.001M EDTA
<u>10X - Loading Buffer (DNA)*</u>		
0.25% Bromophenol blue		
0.25% Xylene cyanol		
20% Ficoll Type 400		
0.1M EDTA, pH 8.0		

Ethidium Bromide Staining

EtBr can be premixed with buffers and agarose for use during electrophoresis. Add to agarose only after temperature has fallen below 55°C. Gels can also be stained after electrophoresis in a soaking tray. Use EtBr at a final concentration of 0.1 µg/ml from a 1mg/ml stock solution. Ethidium Bromide is a powerful mutagen. ALWAYS wear gloves, eye protection and protective clothing. Dispose of solutions in accordance with the safety regulations of your institution.

4.3 References

1. Hames, B.D., Rickwood, D. (ed.) (1990). *Gel Electrophoresis of Nucleic Acids. A Practical Approach*. 2nd edn. IRL Press, Oxford. Ch. 2.
2. Sambrook, J., Fritsch, E.F., Maniatis, T. (1989). *Molecular Cloning. A Laboratory Manual*. 2nd edn. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York. Ch 6.
3. Ausubel, F.M., Brent, R., Kingston, R.E., Moore, D.D., Seidman, J.G., Smith, J.A., Struhl, K. (ed) (1993). *Current Protocols in Molecular Biology*. Vol. 1, Greene Publishing Associates, Inc. and John Wiley & Sons, Inc., Ch. 2.

SECTION 5 Maintenance of Equipment

5.1 Care and Handling

The plastic components of the Midi-Horizontal units are fabricated from acrylic and polycarbonate. Electrodes and connectors are made from pure platinum, stainless steel, and chrome plated brass. As with any laboratory instrument, adequate care ensures consistent and reliable performance.

After each use, rinse buffer chamber, gel tray and combs with de-ionized water. Wipe dry with a soft cloth or paper towel, or allow to air dry. Whenever necessary, all components may be washed gently with water and a non-abrasive detergent, and rinsed and dried as above. *Never* use abrasive cleaners, glass cleaning sprays or scouring pads to clean the components, as these will damage the unit and components.

Additional precautions:

- Do not autoclave or dry-heat sterilize the apparatus or components.
- Do not expose the apparatus or components to phenol, acetone, benzene, halogenated hydrocarbon solvents or alcohols.
- Avoid prolonged exposure of the apparatus or components to UV light.
- Do NOT treat with diethylpyrocarbonate (DEPC)-treated water for extended periods at 37°C. A brief rinse with DEPC-water is sufficient after a thorough wash.

5.2 Maintenance

The following inspection and maintenance procedures will help maintain the safety and reliable performance of the Midi-Horizontal systems. Replacement parts can be ordered by calling

1-858-755-4959 or by contacting your local distributor.

- Banana plugs and power cords should be inspected regularly. If the banana plugs become loose or do not feel friction tight replace the plugs or power cords.
- Should power cord assemblies (connectors, wire or shrouds) show any signs of wear or damage (e.g. cracks, nicks, abrasions, or melted insulation), replace them immediately.
- The platinum wire is secured to the banana jack by compression between a stainless washer and the jack nut. The nut/washer interface should be tight and free of corrosion.

SECTION 6

Horizontal Midi-Gel Kits and Accessories

Cat.#	Item
GCMGU-402T	Horizontal Midi-Gel Productivity Package CE , with tapered baffles. Gel bed dimensions are 11cm x 11cm. Includes two combs, two UVT gel tray, external gel caster and safety cover. Please specify comb cat. # from list below.
GCMGU-502T	Double-Wide Midi-Gel Productivity Package CE , with tapered baffles. Gel bed dimensions are 14cm x 10cm. Includes two combs, two UVT gel tray, external gel caster and safety cover. Please specify comb cat. # from list below.
GCMGU-602T	Horizontal Midi-Gel Productivity Package CE , with tapered baffles. Gel bed dimensions are 14cm x 14cm. Includes two combs, two UVT gel tray, external gel caster and safety cover. Please specify comb cat. # from list below.

Color Coded Midi-Gel Option:

Available for the MGU-402T, MGU-502T and MGU-602T models. To specify color desired, please add the following color code to the end of the catalog number: -Y (Yellow), -P (Purple), -B (Blue), -PK (Pink), or -G (Green). For example, to specify a green MGU-402T place order for catalog # MGU-402T-G.

MIDI-GEL ACCESSORIES

MGU-402T combs (Available in Polycarbonate, or Teflon-Coated Aluminum)**

Cat.#	Comb Dimensions
-------	-----------------

MBU-1005	1mm x 5 well
MBU-1010	1mm x 10 well**
MBU-1012	1mm x 12 well
MBU-1018	1mm x 18 well
MBU-1505	1.5mm x 5 well
MBU-1510	1.5mm x 10 well**
MBU-1512	1.5mm x 12 well
MBU-1518	1.5mm x 18 well

** To specify Teflon combs add "-TF" to end of cat. #.

MGU-502T and MGU-602T combs (Available in Polycarbonate, or Teflon-Coated Aluminum)**

Cat.#	Comb Dimensions
-------	-----------------

SGC14-1001	1mm x 1 well
SGC14-1003E	1mm x 3 equal wells
SGC14-1003M	1mm x 2 well + middle marker
SGC14-1005	1mm x 5 well
SGC14-1008	1mm x 8 well
SGC14-1010	1mm x 10 well
SGC14-1012	1mm x 12 well
SGC14-1013	1mm x 13 well*
SGC14-1016	1mm x 16 well
SGC14-1020	1mm x 20 well
SGC14-1026	1mm x 26 well
SGC14-1027	1mm x 27 well*
SGC14-1501	1.5mm x 1 well
SGC14-1503E	1.5mm x 3 equal wells
SGC14-1503M	1.5mm x 2 well + middle marker
SGC14-1505	1.5mm x 5 well
SGC14-1508	1.5mm x 8 well
SGC14-1510	1.5mm x 10 well
SGC14-1512	1.5mm x 12 well
SGC14-1513	1.5mm x 13 well*
SGC14-1516	1.5mm x 16 well
SGC14-1520	1.5mm x 20 well
SGC14-1526	1.5mm x 26 well
SGC14-1527	1.5mm x 27 well

Cat.#	Comb Dimensions
-------	-----------------

SGC14-2001	2mm x 1 well
SGC14-2003E	2mm x 3 equal wells
SGC14-2003M	2mm x 2 well + middle marker
SGC14-2005	2mm x 5 well
SGC14-2008	2mm x 8 well
SGC14-2010	2mm x 10 well
SGC14-2012	2mm x 12 well
SGC14-2013	2mm x 13 well*
SGC14-2016	2mm x 16 well
SGC14-2020	2mm x 20 well
SGC14-2026	2mm x 26 well
SGC14-2027	2mm x 27 well*
SGC14-3001	3mm x 1 well
SGC14-3003E	3mm x 3 equal wells
SGC14-3003M	3mm x 2 well + middle marker
SGC14-3005	3mm x 5 well
SGC14-3008	3mm x 8 well
SGC14-3010	3mm x 10 well
SGC14-3012	3mm x 12 well
SGC14-3013	3mm x 13 well*
SGC14-3016	3mm x 16 well
SGC14-3020	3mm x 20 well
SGC14-3026	3mm x 26 well
SGC14-3027	3mm x 27 well

*Designed for multi-channel pipettor loading.

** To specify Teflon combs add "-TF" to end of cat. #.

MIDI-GEL ACCESSORIES-CONTINUED

Midi Gel Trays (UV Transparent)

Cat. #	Item
--------	------

MGT-400	Gel Tray for the MGU-402T, 2 slots. Pathlength of 4cm/slot
MGT-500	Gel Tray for the MGU-502T, 2 slots. Pathlength of 3.6cm/slot
MGT-600	Gel Tray for the MGU-602T, 4 slots. Pathlength of 3.0cm/slot

Replacement Baffles

Cat. #	Item
BAFFLES-402T	Tapered Baffles for MGU-402T
BAFFLES-502T	Tapered Baffles for MGU-502T
BAFFLES-602T	Tapered Baffles for MGU-602T

Horizontal Gel Casters

Cat. #	Item
GCT10-111	Gel Casting Tray for 10.5cm (w) x 11cm (l), single place
GCT10-112	Gel Casting Tray for 10.5cm (w) x 11cm (l), two place
GCT14-101	Gel Casting Tray for 14cm (w) x 10cm (l), single place
GCT14-102	Gel Casting Tray for 14cm (w) x 10cm (l), two place
GCT14-161	Gel Casting Tray for 14cm (w) x 16cm (l), single place
GCT14-162	Gel Casting Tray for 14cm (w) x 16cm (l), two place

Gel Tape

Cat. #	Item
GT-72-10	Gel tape, 72 yard x 1.0"
GT-72-15	Gel tape, 72 yard x 1.5"

Power Supply

Cat. #	Item
EPS-200-II	Mini Power Supply, 24 hour timer, CV or CC, 10-200V, 110V/60Hz, current range:4-2000mA, 200 watts
EPS-200-IIV	Mini Power Supply, 24 hour timer, CV or CC, 10-200V, 220V/50Hz, current range:4-2000mA, 200 watts
EPS-300-II	Mini Power Supply, 24 hour timer, CV or CC, 10-300V, 110V/60Hz, current range:4-500mA, 90 watts
EPS-300-IIV	Mini Power Supply, 24 hour timer, CV or CC, 10-300V, 220V/50Hz, current range:4-500mA, 90 watts
EPS-3000-III	Constant Power Supply, 0-3000 Volt, 110V/60Hz, Series II
EPS-3000-IIV	Constant Power Supply, 0-3000 Volt, 220V/50Hz, Series II

Replacement Safety Covers with attached power leads

Cat. #	Item
MGU-402T-SC	Replacement safety cover with power leads for MGU-402T
MGU-502T-SC	Replacement safety cover with power leads for MGU-502T
MGU-602T-SC	Replacement safety cover with power leads for MGU-602T

NOTES

NOTES

CONTACT INFORMATION

Telephone:
Local or International
858-755-4959
Toll Free: 800-243-4959

Fax: 858-755-0733

Online ordering:
www.cbsscientific.com

E-mail address:
sales@cbssci.com

Mailing address:
C.B.S. Scientific Company
P.O. Box 856
Del Mar, CA 92014

Shipping address:
C.B.S. Scientific Company
420 South Cedros
Solana Beach, CA 92075

Credit Card Options:
Visa/Mastercard
Discover/American Express