INSTRUCTION MANUAL

Linear Gradient Makers

GM-20

GM-40

GM-100

GM-150

GM-200

GM-500

GM-1000

GM-2000

TABLE OF CONTENTS

| | Page |
|--|------------|
| Important User Information | 3-4 |
| Section 1 General Information | |
| 1.1 Introduction | 5 |
| 1.2 Specifications | 5 |
| 1.3 Safety | 5 |
| Section 2 Description of parts | |
| 2.1 Unpacking | 6 |
| 2.2 Components/Assembly | 6 |
| Section 3 Instructions for Use | |
| 3.1 Gradient Maker Preparation | 6-7 |
| 3.2 General Notes | 8 |
| Section 4 Maintenance of Equipment | |
| 4.1 Care and Handling | 9 |
| 4.2 Maintenance | 9 |
| Section 5 Accessories | 9 |

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. Scientific warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. Scientific warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. Scientific will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. Scientific shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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
SECTION 1 General Information

1.1 Introduction


The C.B.S. Scientific line of linear gradient makers can be used to make reproducible gradients easily up to 2000mls total volume by filling each cylindrical chamber. Eight different sizes of linear gradient makers with side outlets are available. Gradient formation is controlled by a Teflon valve centered between the two chambers. Gradient flow rate is controlled either by the valve fitted on the outflow male luer adapter, or by a peristaltic pump.

1.2 Specifications

The smaller gradient makers (GM-20/40/100) have a support rod for attachment to a ring stand and come with 10, 20, or 50ml chambers per side (Fig. 1).

The larger sizes are available with 75,100, 250, 500 or 1000mls per chamber arranged side-by-side on a self-supporting base (Fig. 2).

1.3 Safety

Be sure that all gradient makers are secure either by attaching support rod to a ring stand or mounting on a level surface.

SECTION 2 Description of Parts

2.1 Unpacking

Remove the gradient maker from the shipping container and place on a level surface in an authorized area. Verify that the catalog number on the gradient maker matches the packing list. Check for any damage which may have occurred during shipping. Contact C.B.S. Scientific customer service immediately if any damage is found.

2.2 Components/Assembly

The gradient makers consist of the acrylic/teflon mixing chamber made either from a solid block or tubing. Each gradient maker is supplied with an outflow control valve, which should be attached to the male luer adapter on the side of the base.

SECTION 3 Instructions for Use

3.1 Gradient Maker Preparation

1. Prepare gel sandwich and gel solutions as described. See Figure 3 for general set-up of the gel casting procedure.

Vertical Gradient Gel Casting using Gel Wrap, GM-40 gradient former and mini-pump


Figure 3

1. Assemble gel plate sandwich using Gel Wrap®.
2. Stand upright and level on two bottom casting clamps (GPC-0002).
3. Pump desired gradient into plate sandwich at medium speed. (Determine optimal speed empirically.)
4. Insert 16 well comb; allow time for polymerization.
5. Rinse cassette with d.i. H₂O. Remove comb.
6. Transfer sandwich to gel cassette; use GPC-0001 spring clamp. Leave gel wrap in place.

2. When using the 10ml, 20ml or 50ml chambers attach gradient maker to a ring stand support and drop a small magnetic flea or spin-bar into chamber on the right hand side ('B'). Ensure both valves are closed. Set up magnetic stirrer.
3. Outlet # 1 should be fitted with luer type valve and needle with attached tubing to deliver acrylamide mix to gel plate sandwich.
4. Pour solution 'B' into right side of gradient maker. Turn on magnetic stirrer so that it gently mixes the solution.
5. Open inside (middle) valve #2 to allow air bubble to escape. Allow approximately 1ml of solution 'B' to backflow into the left chamber. Remove with pipette and add back to right chamber B.
6. Add solution 'A' to left chamber of gradient maker.
7. Open valve #2 and, then, outside valve #1 to the pump.
8. Turn on peristaltic pump to drive the gel mix into the gel plate sandwich. The flow rate should allow the gel formation to be complete in 10-20 minutes.
9. Exit tube should be attached to near side of gel plate with tape.
10. When all of the gel mix has entered the gel plate sandwich, place the end of outlet tubing from valve #1 into a beaker and flush the chambers and tubing with distilled water so that acrylamide will not harden in the apparatus.

Note: If gel volume is not enough to fill sandwich, use 0% to "top-off" if using the back side of the sawtooth comb for flat interface formation. If using water-saturated Butanol for overlay, leave a 0.5cm void to create a flat interface.

| <u>Gel size</u> | <u>Spacers</u> | <u>Volume required</u> |
|-----------------|-----------------|-----------------------------|
| 16.5cm x 200mm | (yellow) 0.75mm | 20mls (10.5mls per chamber) |
| 16.5cm x 200mm | (red) 1.0mm | 30mls (15.5mls per chamber) |
| 16.5cm x 200mm | (blue) 1.5mm | 40mls (21.0mls per chamber) |

When inserting the flat side of the sawtooth comb, be careful not to exceed 3/16" depth (about half of the tooth length). Otherwise, you will find it difficult to load your sample as there will not be enough "notch" left above the edge of the glass plate. Allow gel to polymerize for 20-30 minutes.

3.2 General Notes

In general, a gradient maker consists of two cylindrical chambers, with connecting channel between them. The two chambers may be arranged side by side or in a concentric arrangement. Each chamber differs according to which flow outlet or “valve” is connected to it. One flow valve is connected between the two chambers. The second flow valve is to be connected to tubing.

The gradient gel should be separated into two solutions; each one-half the final volume of the total gel volume. Each one-half volume should be the high and the low extreme of the gradient. For example, if the gradient is 20 to 50% UF and the total volume of the gel cast is 30mls, then a 15ml solution of 20% UF and 15mls of 50% UF should be prepared.

Make sure that the output tube is secured (tape) to the gel cast sandwich, ready for pouring. If the gel is to be poured by gravity, rather than by pumping, close valve #1.

First, identify which chamber is connected to the flow outlet that connects to the output tubing: this chamber holds the high gradient solution.

Next, ensure that the valve for the flow between the high and low chamber is closed. This valve should not be open until both solutions are poured. If it is opened before both solutions are poured, then the volume poured will have one half of its contents flow into the other chamber, equalizing the volume of solution in each chamber.

Pour one of the solutions into its appropriate chamber.

Flip the valve connecting the two chambers open and close very quickly (without pausing). This serves to get rid of any air bubbles that may be in the valve.

Pour the other solution into its appropriate chamber.

Open the valve connecting the two chambers.

Turn on the pump that controls the output. If the gel is poured through gravity, then open the valve on the output tubing.

If you are pouring a 2D gel, do not let the gel pour all the way to the top of the gel cast; otherwise, there will be no room to load your sample.

If you are doing a 1D gel, it may prove useful to leave 0.75mm of space at the top of the gel and add a neutral Polyacrylamide plug for the comb.

After pouring the gel, use a syringe to put some buffer at the top of the gel. This allows for the gel to settle and polymerize with a flat surface for loading. Disconnect the output tubing from the gel.

If there is excess solution, discard into a beaker through the tubing formerly connected to the gel.

When the gradient maker is empty, add water to clean out the remaining solution and left over polymerized patches.

After polymerization, transfer to gel apparatus for loading.

SECTION 4

Maintenance of Equipment

4.1 Care and Handling

Do not use organic solvents to clean the gradient maker. Use only deionized water. After each use, immediately pump deionized water through both chambers to clear residual acrylamide. Allow to air dry; do not use drying oven.

4.2 Maintenance

Replacement parts can be ordered by calling 1-858-755-4959 or by contacting your local distributor.

SECTION 5

Gradient Makers, Linear

| Cat.# | Item |
|---------|---|
| GM-20 | Gradient Maker, 10ml per side, with side outlet |
| GM-40 | Gradient Maker, 20ml per side, with side outlet |
| GM-100 | Gradient Maker, 50ml per side, with side outlet |
| GM-150 | Gradient Maker, 75ml per side, with side outlet |
| GM-200 | Gradient Maker, 100ml per side, with side outlet |
| GM-500 | Gradient Maker, 250ml per side, with side outlet |
| GM-1000 | Gradient Maker, 500ml per side, with side outlet |
| GM-2000 | Gradient Maker, 1000ml per side, with side outlet |

NOTES

NOTES

CONTACT US

Telephone:
Local or International
858-755-4959
Toll Free: 858-755-0733

Fax: 858-755-0733

Check out our new and improved online ordering!

www.cbsscientific.com

E-mail address:
sales@cbssci.com

Mailing address:
C.B.S. Scientific Company
P.O. Box 856
Del Mar, CA 92014

Shipping address:
C.B.S. Scientific Company
10805 Vista Sorrento Pkwy
Suite 100
San Diego, CA 92121

Credit Card Options:
Visa/Mastercard
Discover/American Express